

Chloral Hydrate Solution (10% w/v)

Cat #: A-CSH489

Size: 100mL, 500mL

Storage: 2~8°C, avoid light (6 months)

Product Description

The Chloral hydrate solution (10% w/v) is composed of 10% Chloral hydrate and deionized water. It has an irritating odor and corrosive properties, so it should be stored in a sealed manner. This reagent is for research purposes only and should not be used for clinical diagnosis or other purposes.

Usage

Follow the specific experimental requirements.

Notes

1. Store in a sealed manner to prevent volatilization and corrosion of other substances.
2. For your safety and health, please wear lab coat and disposable gloves during operation.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.