

PH Meter Electrode Filling Solution

Cat #: A-CSH492

Size: 100mL

Storage: Room temperature (36 months)

Product Description

Prolonged exposure to acid and alkali or improper storage solutions can cause aging of the pH electrode. The pH Electrode Refilling Solution effectively protects the pH electrode.

Usage

Follow the specific experimental requirements. After cleaning the electrode with ultrapure water, store the electrode in the pH Electrode Refilling Solution.

Notes

For your safety and health, please wear lab coat and disposable gloves during operation.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.