

L-Glutamine Solution (200mmol/L, Sterile Filtered)

Cat #: A-CSH493

Size: 100mL

Storage: -20°C (12 months)

Product Description

L-Glutamine Solution (200mmol/L) is composed of L-glutamine, sodium chloride, deionized water, and undergoes sterile filtration. It can be used as an additive in mammalian cell culture media and has a wide range of applications.

Usage

Operate according to the specific requirements of the experiment

Notes

1. Avoid repeated freeze-thaw cycles and pay attention to aseptic techniques.
2. L-Glutamine is unstable in solution and can decompose when stored at 4°C for 1 week. Therefore, it should be added to the culture medium before use.
3. If a culture medium containing L-glutamine is stored in a refrigerator at 4°C for more than 2 weeks, the original L-glutamine should be replenished.
4. For your safety and health, please wear laboratory attire and disposable gloves when handling

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.