

Sodium Chloride Solution (5mol/L)

Cat #: A-CSH494

Size: 100mL

Storage: Room temperature (12 months)

Product Description

Sodium chloride is a fundamental component of various research reagents and has extremely wide applications.

Sodium Chloride Solution (5mol/L) is composed of sodium chloride and deionized water.

Usage

Operate according to the specific requirements of the experiment.

Notes

1. Pay attention to aseptic techniques to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves when handling.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.