

Glycerol Solution (100%, Sterile)

Cat #: A-CSH496

Size: 100mL

Storage: Room temperature (12 months)

Product Description

Glycerol Solution (100%, sterile) is a solution containing 100% high-quality glycerol. It has undergone high-pressure sterilization and is commonly used as a solvent, lubricant, and protective agent. This reagent is intended for research purposes only and should not be used for clinical diagnosis or other purposes.

Usage

Dilute to the appropriate concentration before use

Notes

1. Pay attention to aseptic techniques to avoid bacterial contamination.
2. Glycerol Solution (100%, sterile) is quite viscous. When diluting, it is recommended to first dilute it significantly to obtain a dilute solution before use. Otherwise, there may be differences in concentration.
3. For your safety and health, please wear laboratory attire and disposable gloves when handling.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.