

Glycerol Solution (10%, Sterile)

Cat #: A-CSH499

Size: 500mL

Storage: Room temperature (6 months)

Product Description

Glycerol Solution (10%, sterile) is a solution containing 10% high-quality glycerol. It has undergone high-pressure sterilization and is commonly used as a solvent, lubricant, and protective agent. This reagent is intended for research purposes only and should not be used for clinical diagnosis or other purposes.

Usage

Dilute to the appropriate concentration before use

Notes

1. Pay attention to aseptic techniques to avoid bacterial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves when handling.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.