

Ammonium Chloride Solution (50mmol/L)

Cat #: A-CSH504

Size: 100mL

Storage: 2-8°C (6 months)

Product Description

Ammonium Chloride, with a molecular weight of 53.49 and a chemical formula of NH_4Cl , is a strong electrolyte. When dissolved in water, it ionizes into ammonium ions and chloride ions. It can neutralize acidic intracellular components and inhibit the synthesis of nerve sheath amino acids. Ammonium Chloride has multiple applications.

Usage

Follow the specific experimental requirements for operation. The working concentration is usually in the range of 1-50mmol/L.

Notes

For your safety and health, please wear a lab coat and disposable gloves while conducting the procedure.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.