

## Ammonium Chloride Solution (5mol/L)

Cat #: A-CSH505

Size: 100mL

Storage: 2-8°C (6 months)

### Product Description

Ammonium Chloride has a molecular weight of 53.49 and a chemical formula of  $\text{NH}_4\text{Cl}$ . It is a strong electrolyte that dissolves in water, ionizing into ammonium ions and chloride ions. It can neutralize acidic intracellular components and inhibit the synthesis of nerve sheath amino acids. Ammonium Chloride has various applications.

### Usage

Follow the specific experimental requirements for operation. The working concentration is commonly in the range of 1-50mmol/L.

### Notes

For your safety and health, please wear a lab coat and disposable gloves while performing the procedure.

### Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.