

Glycerol Solution (50%, Sterile)

Cat #: A-CSH510

Size: 100mL

Storage: Room temperature (12 months)

Product Description

Glycerol Solution (50%, Sterile) is a solution of high-quality glycerol at a concentration of 50%. It has been subjected to high-pressure sterilization and is commonly used as a solvent, lubricant, and protective agent. This reagent is for research use only and is not intended for clinical diagnosis or other purposes.

Usage

Dilute to the appropriate concentration before use.

Notes

1. Pay attention to aseptic operation to avoid contamination by bacteria or other microorganisms.
2. For your safety and health, please wear a lab coat and disposable gloves while performing the procedure.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.