

PIPES Buffer (1mol/L, pH 6.8)

Cat #: A-CSH519

Size: 100mL

Storage: 2-8°C (12 months)

Product Description

PIPES is a non-ionic zwitterionic buffer with good buffering capacity in the pH range of 6.1 to 7.5. PIPES Solution (1mmol/L, pH 6.8) is mainly composed of PIPES and deionized water. It is a commonly used biological buffer solution and is often used in flow cytometry to analyze the cell cycle.

Usage

Operate according to specific experimental requirements.

Notes

1. High concentrations of PIPES may have some toxicity to certain cells.
2. For your safety and health, please wear lab coat and disposable gloves during operation.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.