

Bovine Serum Albumin Solution (1% BSA)

Cat #: A-CSH528

Size: 100mL

Storage: 2-8°C (6 months)

Product Description

Bovine Serum Albumin Solution (1% BSA) is a reagent used for coating and protein supplementation, among other purposes. It mainly consists of BSA and deionized water.

Usage

Operate according to the specific requirements of the experiment.

Notes

1. Handle with care to prevent microbial contamination.
2. For your safety and health, please wear a lab coat and disposable gloves while handling.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.