

Sodium Phosphate Buffer (0.01mol/L, pH 5.8-8.0)

Cat #: A-CSH530

Size: 500mL

Storage: 2~8°C (12 months)

Product Description

Sodium phosphate buffer (0.01mol/L, pH 5.8-8.0) is prepared by mixing disodium hydrogen phosphate and monosodium phosphate in different proportions. The phosphate ion concentration in the buffer is 0.01mol/L, and the customer can choose the pH value according to their specific needs.

Usage

Operate according to the specific requirements of the experiment.

Notes

For your safety and health, please wear a lab coat and disposable gloves while handling.

Disclaimer

For your safety and health, please wear a lab coat and disposable gloves while handling.