

Sodium Phosphate Buffer (1mol/L, pH 5.8-8.0)

Cat #: A-CSH533

Size: 500mL

Storage: 2~8°C (12 months)

Product Description

Sodium Phosphate Buffer (1mol/L, pH 5.8-8.0) is prepared by mixing monosodium phosphate and disodium phosphate in different proportions. The concentration of phosphate ions is 1mol/L, and customers can choose the pH value within the range of 5.8 to 8.0 according to their needs.

Usage

Follow the specific requirements of your experiment.

Notes

For your safety and health, please wear lab attire and disposable gloves while handling.

Disclaimer

For your safety and health, please wear lab attire and disposable gloves while handling.