

Sodium Phosphate Buffer (0.05mol/L, pH 7.8)

Cat #: A-CSH535

Size: 500mL

Storage: 2~8°C (12 months)

Product Description

Sodium Phosphate Buffer (0.05mol/L, pH 7.8) is prepared by mixing monosodium phosphate and disodium phosphate in different proportions. The concentration of phosphate ions is 0.05mol/L, and customers can choose the pH value according to their needs.

Usage

Follow the specific requirements of your experiment.

Notes

For your safety and health, please wear lab attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.