

Sodium Phosphate Buffer (50mmol/L, with EDTA)

Cat #: A-CSH537

Size: 500mL

Storage: 2~8°C (12 months)

Product Description

Phosphate Buffer (50mmol/L, containing EDTA) is composed of 50mM Na₂HPO₄, 300mM NaCl, and 2mM EDTA. It has a pH value of 6.0~6.5. This buffer can be used for protein purification experiments. Please note that this reagent is for research purposes only and is not suitable for clinical diagnosis or other applications.

Usage

Dilute to the appropriate concentration according to the specific experimental requirements before use.

Notes

For your safety and health, please wear lab attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.