

Borate Buffer (0.002M, pH 8.0)

Cat #: A-CSH550

Size: 500mL

Storage: Room temperature (12 months)

Product Description

Borate buffer is composed of boric acid and sodium chloride, with a boric acid concentration of 0.002mol/L. The pH is adjusted to 8.0 using NaOH.

Usage

Follow the specific requirements of your experiment.

Notes

1. It may cause irritation to the human body, so please take appropriate precautions.
2. For your safety and health, please wear lab attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.