

Na₂HPO₄-KH₂PO₄ Buffer (0.2mol/L, pH 5.8-8.0)

Cat #: A-CSH551

Size: 500mL

Storage: 2~8°C (12 months)

Product Description

Sodium Dihydrogen Phosphate-Potassium Dihydrogen Phosphate Buffer (0.2mol/L, pH 5.8~8.0) is prepared by mixing Sodium Dihydrogen Phosphate and Potassium Dihydrogen Phosphate in different proportions. The concentration of phosphate ions is 0.2mol/L. Customers can choose from pH values ranging from 5.8 to 8.0 according to their specific needs.

Usage

Follow the specific requirements of your experiment.

Notes

1. It may cause irritation to the human body, so please take appropriate precautions.
2. For your safety and health, please wear lab attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.