

Na₂HPO₄-KH₂PO₄ Buffer (25mmol/L, pH 7.3)

Cat #: A-CSH554

Size: 500mL

Storage: 2~8°C (12 months)

Product Description

Sodium dihydrogen phosphate-potassium dihydrogen phosphate buffer (25mmol/L, pH 7.3) is a mixture of sodium dihydrogen phosphate and potassium dihydrogen phosphate in different proportions. The concentration of phosphate ions is 0.2mol/L. Customers can choose the pH value according to their specific needs, ranging from pH 5.8, 5.9, 6.0, 6.1, 6.2, 6.3, 6.4, 6.5, 6.6, 6.7, 6.8, 6.9, 7.0, 7.1, 7.2, 7.3, 7.4, 7.5, 7.6, 7.7, 7.8, 7.9, to 8.0.

Usage

Operate according to specific experimental requirements.

Notes

For your safety and health, please wear laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.