

Acetic Acid Solution (0.5%)

Cat #: A-CSH563

Size: 500mL

Storage: Room temperature (12 months)

Product Description

Acetic acid aqueous solution (0.5%) is primarily composed of glacial acetic acid and deionized water. It is an essential auxiliary reagent.

Usage

Operate according to specific experimental requirements.

Notes

1. Store in a sealed container and use it as soon as possible after opening, as the active component is volatile.
2. For your safety and health, wear laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.