

## Acetic Acid Solution (1%)

Cat #: A-CSH564

Size: 500mL

Storage: Room temperature (12 months)

### Product Description

The 1% acetic acid aqueous solution is primarily composed of glacial acetic acid and deionized water. It is an essential auxiliary reagent.

### Usage

Operate according to specific experimental requirements.

### Notes

1. Store in a tightly sealed container and use it as soon as possible after opening, as the active components are volatile.
2. For your safety and health, wear laboratory attire and disposable gloves when handling.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.