

Sodium Acetate Buffer (0.1mol/L, pH 6.1)

Cat #: A-CSH573

Size: 500mL

Storage: Room temperature (12 months)

Product Description

Sodium acetate solution, also known as sodium ethanoate solution, is primarily composed of NaAc and glacial acetic acid. It mainly provides acetate ions and helps regulate pH value.

Usage

Use according to specific experimental requirements.

Notes

1. Dispense according to the required amount.
2. Store in a tightly sealed container to prevent an increase in pH value.
3. For your safety and health, wear laboratory attire and disposable gloves when handling..

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.