

## NaCO<sub>3</sub>-NaHCO<sub>3</sub> Buffer (0.1mol/L, pH 9.16-10.83)

Cat #: A-CSH577

Size: 500mL

Storage: Room temperature (12 months)

### Product Description

Sodium Carbonate-Sodium Bicarbonate Buffer Solution (0.1mol/L, pH 9.16-10.83) is composed of sodium carbonate and sodium bicarbonate. Customers can choose from the following pH values: 9.16, 9.4, 9.51, 9.79, 9.9, 10.14, 10.28, 10.53, 10.83. It is important to store the solution in a sealed container.

### Usage

Use according to specific experimental requirements.

### Notes

1. Store the solution in a sealed container to prevent evaporation.
2. For your safety and health, please wear a lab coat and disposable gloves when handling.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.