

Glucose Solution (50%, Sterile)

Cat #: A-CSH578

Size: 100mL, 500mL

Storage: 2-8°C (6 months)

Product Description

Glucose Solution (50%, Sterile) is composed of glucose and deionized water. It has undergone sterile processing and can be used as an additive in cell culture media. It has a wide range of applications and is an essential auxiliary reagent.

Usage

Operate according to specific experimental requirements.

Notes

1. Pay attention to sterile procedures to avoid microbial contamination.
2. For your safety and health, please wear a lab coat and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.