

Standard Glucose Solution (100mmol/L)

Cat #: A-CSH584

Size: 1mL

Storage: 2-8°C (6 months)

Product Description

Standard Glucose Solution (100mmol/L) is composed of glucose and high-purity deionized water. It can be used as a glucose standard.

Usage

1. Operate according to specific experimental requirements.
2. Generally, prepare multiple dilutions to create a standard curve.

Notes

1. Avoid repeated freezing and thawing.
2. For your safety and health, please wear a lab coat and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.