

Sucrose Solution (40%)

Cat #: A-CSH589

Size: 100mL

Storage: 2-8°C (4 months)

Product Description

Sucrose Solution (40%) is mainly composed of sucrose and deionized water. It can be used for tissue dehydration or other purposes and is an essential auxiliary reagent.

Usage

Operate according to specific experimental requirements.

Notes

1. Store in a sealed container and avoid microbial contamination.
2. For your safety and health, please wear a lab coat and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.