

Sucrose Solution (2mol/L)

Cat #: A-CSH592

Size: 500mL

Storage: 2-8°C (4 months)

Product Description

Sucrose solution (2mol/L) is primarily composed of sucrose and deionized water. It can be used for tissue dehydration or other purposes and is an essential auxiliary reagent.

Usage

Operate according to specific experimental requirements.

Notes

1. Store in a sealed container and avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.