

Sodium Acetate Buffer (0.1mol/L, pH 3.6-5.6, Sterile)

Cat #: A-CSH601

Size: 500mL

Storage: Room temperature (12 months)

Product Description

Sodium acetate buffer solution (0.1mol/L, pH 3.6~5.6, sterile) is prepared by mixing acetic acid and sodium acetate in different proportions. The concentration of acetate ions is 0.05mol/L, and it is primarily used as an analytical reagent for pH adjustment. Customers can choose from pH 3.6, 3.8, 4.0, 4.2, 4.4, 4.6, 4.8, 5.0, 5.2, 5.6, or 5.8, and the solution is subjected to high-pressure sterilization.

Usage

Operate according to specific experimental requirements.

Notes

For your safety and health, please wear lab coat and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.