

Acetate Buffer (0.05mol/L, pH 5.0)

Cat #: A-CSH603

Size: 500mL

Storage: Room temperature (12 months)

Product Description

Acetate buffer solution (0.05mol/L, pH 5.0) is composed of sodium acetate and glacial acetic acid. It primarily provides acetate ions and is used as an analytical reagent for pH adjustment.

Usage

Operate according to specific experimental requirements.

Notes

For your safety and health, please wear lab coat and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.