

Sodium Acetate Solution (3mol/L, pH 5.0, Sterile)

Cat #: A-CSH609

Size: 500mL

Storage: Room temperature (12 months)

Product Description

NaAc (3mol/L, pH 5.0, sterile) is mainly composed of 3M NaAc, adjusted to pH 5.0. It is commonly used for DNA ethanol precipitation, etc.

Usage

Operate according to specific experimental requirements.

Notes

1. If the amount used each time is small, it can be appropriately aliquoted before use.
2. For your safety and health, please wear lab coat and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.