

Sodium Acetate Solution (3mol/L, pH 5.2, Sterile)

Cat #: A-CSH610

Size: 100mL, 500mL

Storage: Room temperature (12 months)

Product Description

Sodium Acetate Solution (3M, pH 5.2, sterile) is primarily composed of 3M NaAc, adjusted to pH 5.2, and is commonly used for DNA ethanol precipitation, among other applications.

Usage

Operate according to specific experimental requirements.

Notes

1. If the required amount for each use is small, it can be aliquoted appropriately before use.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.