

EDTA Solution (1mol/L, pH 7.4)

Cat #: A-CSH623

Size: 500mL

Storage: Room temperature (12 months)

Product Description

EDTA solution (1 M, pH 7.4) is prepared using EDTA•2Na and adjusted to pH 7.4. This solution is commonly used for chelating metal ions and is a frequently used reagent in molecular biology.

Usage

Operate according to specific experimental requirements.

Notes

1. When sterile conditions are required for the experiment, please pay attention to aseptic techniques.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.