

Calcium Chloride Solution, CaCl_2 (0.005mol/L, Sterile)

Cat #: A-CSH647

Size: 100mL

Storage: Room temperature (6 months)

Product Description

Calcium chloride is a fundamental component of various research reagents and can also be used as a desiccant. The calcium chloride solution (0.005mol/L, sterile) is composed of calcium chloride and deionized water, and it has undergone high-pressure sterilization.

Usage

Operate according to specific experimental requirements.

Notes

1. Pay attention to sterile operations to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.