

Calcium Chloride Solution, CaCl_2 (1mol/L, Sterile)

Cat #: A-CSH650

Size: 100mL, 500mL

Storage: Room temperature (6 months)

Product Description

Calcium chloride is a fundamental component of various research reagents and can also be used as a drying agent, among other applications. The Calcium Chloride Solution (1mol/L, Sterile) is composed of calcium chloride and deionized water, and it undergoes high-pressure sterilization.

Usage

Operate according to the specific requirements of the experiment.

Notes

1. Pay attention to aseptic techniques to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves during operation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.