

Sodium Chloride Solution, NaCl (0.1mol/L, Sterile)

Cat #: A-CSH662

Size: 100mL

Storage: Room temperature (12 months)

Product Description

Sodium chloride is a fundamental component of various research reagents and has a wide range of applications.

Sodium chloride solution (0.1 mol/L, sterile) is composed of sodium chloride and deionized water, and it undergoes high-pressure sterilization.

Usage

Operate according to specific experimental requirements.

Notes

1. Pay attention to aseptic techniques to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves during the operation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.