

Sodium Chloride Solution, NaCl (1mol/L, Sterile)

Cat #: A-CSH668

Size: 100mL, 500mL

Storage: Room temperature (12 months)

Product Description

Sodium chloride is a fundamental component of various research reagents and has a wide range of applications. The sodium chloride solution (1 M, sterile) is composed of sodium chloride and deionized water, and it has undergone high-pressure sterilization treatment.

Usage

Operate according to the specific requirements of the experiment.

Notes

1. Pay attention to aseptic operation to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.