

## Sodium Chloride Solution, NaCl (3mol/L, Sterile)

Cat #: A-CSH670

Size: 100mL

Storage: Room temperature (12 months)

### Product Description

Sodium chloride is a fundamental component of various research reagents and has a wide range of applications. The sodium chloride solution (3 M, sterile) is composed of sodium chloride and deionized water, and it has undergone high-pressure sterilization treatment.

### Usage

Operate according to the specific requirements of the experiment.

### Notes

1. Pay attention to aseptic operation to avoid microbial contamination.
2. For your safety and health, please wear laboratory attire and disposable gloves while handling.

### Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.