

Calcium Chloride Solution, CaCl_2 (6.8mol/L, Sterile)

Cat #: A-CSH732

Size: 100mL

Storage: Room temperature (6 months)

Product Description

Calcium chloride is a basic component of various research reagents and can also be used as a drying agent. The calcium chloride solution (6.8mol/L, sterile) is composed of calcium chloride and deionized water, and has been subjected to high-pressure sterilization.

Usage

Operate according to the specific requirements of the experiment.

Notes

1. Pay attention to aseptic operation to avoid microbial contamination.
2. For your safety and health, please wear lab coat and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.