

Sterile Ultra-Pure Water

Cat #: A-CSH749

Size: 500mL

Storage: Room temperature (12 months)

Product Description

Sterile Ultrapure Water, also known as UP water, refers to water with a resistivity of $18 \text{ M}\Omega \cdot \text{cm}$ (25°C) and minimal impurities. It is primarily used for the preparation of various culture media and solutions in cell culture experiments.

Usage

Use according to the specific experimental requirements.

Notes

1. Take precautions to prevent microbial contamination of the solution.
2. This product is intended for scientific research purposes by professionals only and should not be used for clinical diagnosis or treatment, or for food or drug purposes.
3. For your safety and health, please wear laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.