

DMEM/F12 Medium (With Penicillin-Streptomycin)

Cat #: A-CSH767

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM/F-12 (Dulbecco's Modified Eagle Medium/Nutrient Mixture F-12) is a widely used basal medium for supporting the growth of various mammalian cells, including MDCK, glial cells, fibroblasts, human endothelial cells, and rat fibroblasts. In addition to the incomplete DMEM/F-12 medium, we also offer DMEM/F-12 medium supplemented with antibiotics and fetal bovine serum to meet different customer needs.

This medium contains: 3.15 g/L glucose, 0.365 g/L L-glutamine, phenol red, sodium bicarbonate, 50KU penicillin, 50 mg streptomycin, etc. This medium contains 3.57 g/L HEPES.

Usage

Use according to specific experimental requirements and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination of the solution and avoid freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and should not be used for clinical diagnosis or treatment, or for food or drugs.
3. For your safety and health, please wear lab attire and disposable gloves during operation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.