

DMEM Medium (High Glucose, With FBS, Penicillin-Streptomycin)

Cat #: A-CSH772

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM (Dulbecco's Modified Eagle Medium) is a widely used basal medium for supporting the growth of various mammalian cells, including primary fibroblasts, neurons, glial cells, HUVECs, smooth muscle cells, as well as cell lines such as HeLa, 293, COS-7, and PC-12. In addition to the incomplete DMEM medium (high glucose), we also provide DMEM medium supplemented with antibiotics and fetal bovine serum to meet different customer needs.

This medium contains: 4.5g/L glucose, 0.584g/L L-glutamine, phenol red, sodium pyruvate, sodium bicarbonate, 10% fetal bovine serum, 50KU penicillin, 50mg streptomycin, etc. This medium does not contain HEPES.

Usage

Use according to specific experimental requirements and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination of the solution and avoid freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and should not be used for clinical diagnosis or treatment, or for food or drugs.
3. For your safety and health, please wear lab coat and disposable gloves during operation.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.