

Ham'S F-12K Medium (With FBS, Penicillin-Streptomycin)

Cat #: A-CSH775

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

Ham's F-12K (Ham's F-12K (Kaighn's) Medium) is an improved version of Ham's F-12K, with increased concentrations of amino acids and pyruvate, reduced glucose dosage, and modified salt composition and content compared to Ham's F-12. Ham's F-12K medium was originally designed for the culture of primary human hepatocytes and differentiated cells of rats and chickens. In addition to the incomplete Ham's F-12K medium, we also provide Ham's F-12K medium supplemented with fetal bovine serum and dual antibodies to meet different customer needs.

This medium contains 1.26 g/L glucose, 0.292g/L L-glutamine, phenol red, sodium bicarbonate, 10% fetal bovine serum, 50KU penicillin, and 50mg streptomycin. This medium does not contain HEPES.

Usage

Use according to specific experimental requirements, and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination and freeze-thaw cycles.
2. This product is intended for scientific research by professionals only and should not be used for clinical diagnosis or treatment, or for food or drugs.
3. For your safety and health, please wear appropriate laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.