

DMEM Medium (Low Glucose, With Penicillin-Streptomycin)

Cat #: A-CSH779

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM (Dulbecco's Modified Eagle Medium) is a widely used basal medium that supports the growth of various mammalian cells, including primary fibroblasts, neurons, glial cells, HUVECs, smooth muscle cells, and cell lines such as HeLa, 293, COS-7, and PC-12. In addition to the incomplete DMEM medium (low glucose), we also provide DMEM medium supplemented with dual antibodies and fetal bovine serum to meet different customer needs.

This medium contains 1.0g/L glucose, 0.584g/L L-glutamine, phenol red, sodium pyruvate, sodium bicarbonate, 50KU penicillin, and 50mg streptomycin. Take precautions to avoid microbial contamination and freeze-thaw cycles.

Usage

Use according to specific experimental requirements, and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination and freeze-thaw cycles.
2. This product is intended for scientific research by professionals only and should not be used for clinical diagnosis or treatment, or for food or drugs.
3. For your safety and health, please wear appropriate laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.