

RPMI 1640 Medium

Cat #: A-CSH783

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

RPMI 1640 medium was originally developed for the suspension culture of monolayer human leukemia cells. Roswell Park Memorial Institute (RPMI) 1640 medium has since been widely used for the cultivation of various mammalian cells, including HeLa, Jurkat, MCF-7, PC12, PBMC, astrocytes, and others. In addition to the regular RPMI 1640 medium, we also offer RPMI 1640 medium supplemented with antibiotics and fetal bovine serum to meet different customer needs. This medium contains 2g/L glucose, 0.3g/L L-glutamine, phenol red, sodium bicarbonate, etc. This medium does not contain HEPES.

Notes

1. When using this product for specific experiments, it is important to take precautions to avoid microbial contamination, and avoid freeze-thaw cycles.
2. This product is intended for scientific research purposes by professionals only and is not to be used for clinical diagnosis or treatment, or for food or drug purposes.
3. For your safety and well-being, please wear laboratory attire and disposable gloves when handling this product.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.