

RPMI 1640 Medium (Without Phenol Red)

Cat #: A-CSH789

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

RPMI 1640 medium was originally developed for the suspension culture of human leukemia cells. Roswell Park Memorial Institute (RPMI) 1640 medium has since been widely used for the culture of various mammalian cells, including HeLa, Jurkat, MCF-7, PC12, PBMCs, astrocytes, and more. In addition to the incomplete RPMI 1640 (Phenol Red-Free) medium, we also offer RPMI 1640 medium supplemented with antibiotics and fetal bovine serum to meet different customer needs.

This medium contains 2g/L glucose, 0.3g/L L-glutamine, sodium bicarbonate, and other components. This medium does not contain HEPES or phenol red.

Usage

Use according to specific experimental requirements and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination of the solution and freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and should not be used for clinical diagnosis or treatment, or for food or drug purposes.
3. For your safety and health, please wear laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.