

DMEM Medium (No Glucose, With Phenol Red)

Cat #: A-CSH802

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM (Dulbecco's Modified Eagle Medium) is a widely used basal medium designed to support the growth of various mammalian cells, including primary fibroblasts, neurons, glial cells, HUVECs, smooth muscle cells, as well as cell lines such as HELA, 293, COS-7, and PC-12. This medium contains 0.584g/L L-glutamine, phenol red, sodium pyruvate, sodium bicarbonate, etc. Precautions should be taken to avoid microbial contamination and freeze-thaw cycles.

Usage

Use according to specific experimental requirements, and take precautions to avoid microbial contamination

Notes

1. Take precautions to avoid microbial contamination and freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and is not for clinical diagnosis or treatment, or for use in food or drugs.
3. For your safety and health, please wear appropriate laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.