

RPMI 1640 Medium (No Glucose, With Phenol Red)

Cat #: A-CSH805

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

RPMI 1640 Medium was initially developed for the suspension culture of human leukemia cells. Roswell Park Memorial Institute (RPMI) 1640 Medium has since been widely used for the cultivation of various mammalian cells, including Hela, Jurkat, MCF-7, PC12, PBMC, and astrocytes. In addition to RPMI 1640 Medium (Sugar-Free, Phenol Red-Containing), we also offer RPMI 1640 medium supplemented with antibiotics and fetal bovine serum to meet different customer needs. This medium contains phenol red, 0.3g/L L-glutamine, sodium bicarbonate, etc. It does not contain HEPES or glucose.

Usage

Use according to specific experimental requirements, and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination and freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and is not for clinical diagnosis or treatment, or for use in food or drugs.
3. For your safety and health, please wear appropriate laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.