

DMEM Medium (High Glucose, With Phenol Red , Without Sodium Pyruvate, L-Glutamine)

Cat #: A-CSH806

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM (Dulbecco's Modified Eagle Medium) is a widely used basal medium designed to support the growth of various mammalian cells, including primary fibroblasts, neurons, glial cells, HUVECs, smooth muscle cells, as well as cell lines such as HELA, 293, COS-7, and PC-12. In addition to the DMEM (High Glucose) basal medium, we also offer DMEM (High Glucose) medium supplemented with antibiotics and fetal bovine serum to meet different customer needs.

This medium contains 4.5g/L glucose, phenol red, sodium bicarbonate, etc. It does not contain HEPES, sodium pyruvate, or L-glutamine.

Usage

Use according to specific experimental requirements, and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination and freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and is not for clinical diagnosis or treatment, or for use in food or drugs.
3. For your safety and health, please wear appropriate laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.