

DMEM/F12 Medium (Without Phenol Red)

Cat #: A-CSH808

Size: 500mL

Storage: 2-8°C (12 months)

Product Description

DMEM/F-12 (Dulbecco's Modified Eagle Medium/Nutrient Mixture F-12) is a widely used basal medium designed to support the growth of various mammalian cells, including MDCK, glial cells, fibroblasts, human endothelial cells, and rat fibroblasts. In addition to the incomplete DMEM/F12 medium, we also offer DMEM/F12 medium supplemented with antibiotics and fetal bovine serum to meet different customer needs.

This medium contains 3.15 g/L glucose, 0.365g/L L-glutamine, sodium bicarbonate, etc. It also contains 3.57g/L HEPES.

This medium does not contain phenol red.

Usage

Use according to specific experimental requirements, and take precautions to avoid microbial contamination.

Notes

1. Take precautions to avoid microbial contamination and freeze-thaw cycles.
2. This product is intended for scientific research use by professionals only and is not for clinical diagnosis or treatment, or for use in food or drugs.
3. For your safety and health, please wear appropriate laboratory attire and disposable gloves when handling.

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.