

Rituximab Biosimilar, Research Grade

Summary

Catalog No.	BGT-ABS-197
Alternative Names	IDEC-C2B8, CAS: 174722-31-7
Host species	Chimeric
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG1-kappa
Applications	Research Grade Biosimilar
Target	B-lymphocyte surface antigen B1, Membrane-spanning 4-domains subfamily A member 1, MS4A1, Leukocyte surface antigen Leu-16, B-lymphocyte antigen CD20, Bp35, CD20
Purification	Protein A or G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Expression system	XtenCHO
Accession	P11836


Stability and Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.

Note

For research use only. Not suitable for clinical or therapeutic use.

Data Image



SDS PAGE for Rituximab

SDS-PAGE