

Evolocumab Biosimilar, Research Grade

Summary

Catalog No.	BGT-ABS-342
Alternative Names	AMG 145, Anti-PCSK9, CAS: 1256937-27-5
Clone ID	Evolocumab
Host species	Human
Species reactivity	Human
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2-lambda
Applications	Research Grade Biosimilar
Target	Proprotein convertase subtilisin/kexin type 9, Subtilisin/kexin-like protease PC9, NARC-1, PC9, Proprotein convertase 9, Neural apoptosis-regulated convertase 1, NARC1, PCSK9
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Expression system	Mammalian Cells
Accession	Q8NBP7


Stability and Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Note


For research use only. Not suitable for clinical or therapeutic use.

Data Image



SDS PAGE for Evolocumab

SDS-PAGE



The purity of this product is >95% as determined by SEC-HPLC.

SEC-HPLC